

## Dosiertabelle

| Konz. in %   | 0,25%   | 0,5%   | 0,75%   | 1%     | 1,5%   | 2%     | 2,5%   | 3%     | 4%      | 5%      |
|--------------|---------|--------|---------|--------|--------|--------|--------|--------|---------|---------|
| Liter-Lösung |         |        |         |        |        |        |        |        |         |         |
| 1 Liter      | 2,5 ml  | 5 ml   | 7,5 ml  | 10 ml  | 15 ml  | 20 ml  | 25 ml  | 30 ml  | 40 ml   | 50 ml   |
| 2 Liter      | 5 ml    | 10 ml  | 15 ml   | 20 ml  | 30 ml  | 40 ml  | 50 ml  | 60 ml  | 80 ml   | 100 ml  |
| 3 Liter      | 7,5 ml  | 15 ml  | 22,5 ml | 30 ml  | 45 ml  | 60 ml  | 75 ml  | 90 ml  | 120 ml  | 150 ml  |
| 4 Liter      | 10 ml   | 20 ml  | 30 ml   | 40 ml  | 60 ml  | 80 ml  | 100 ml | 120 ml | 160 ml  | 200 ml  |
| 5 Liter      | 12,5 ml | 25 ml  | 37,5 ml | 50 ml  | 75 ml  | 100 ml | 125 ml | 150 ml | 200 ml  | 250 ml  |
| 8 Liter      | 20 ml   | 40 ml  | 60 ml   | 80 ml  | 120 ml | 160 ml | 200 ml | 240 ml | 320 ml  | 400 ml  |
| 10 Liter     | 25 ml   | 50 ml  | 75 ml   | 100 ml | 150 ml | 200 ml | 250 ml | 300 ml | 400 ml  | 500 ml  |
| 20 Liter     | 50 ml   | 100 ml | 150 ml  | 200 ml | 300 ml | 400 ml | 500 ml | 600 ml | 800 ml  | 1000 ml |
| 30 Liter     | 75 ml   | 150 ml | 225 ml  | 300 ml | 450 ml | 600 ml | 750 ml | 900 ml | 1200 ml | 1500 ml |

Anwendungsbeispiel: 10 Liter : 1 %igen Anwendungslösung = 100 ml  
Desinfektionsmittel + 9,9 Liter Wasser